


## Automatic With Power Locking

Compact, lightweight and automatic, the AL570 is also available as AL570XL with an extra large platform.

### Features & Benefits

- Dock 'N' Lock automatically secures power chair to platform with positive lock device
- Adjusts to accommodate a wide range of power chairs; no chair modifications required
- Manual crank backup and license plate mount included
- Swing-Away option available



Specifications	AL570	AL570XL Extra Large
Accommodates	Power Chairs	Larger Power Chairs
Platform	26.5" x 34.5" / 35" x 34.5"	26.5" x 34.5" / 35" x 34.5"
Lifting Capacity	350 lbs. (160 kg)	350 lbs. (160 kg)
Installed Weight	69 lbs. (31 kg)	79 lbs. (36 kg)
Securement	Dock 'N' Lock	Dock 'N' Lock
Hitch	Class III or greater	Class III or greater
Options	Battery pack, Swing-Away, Power chair covers	Battery pack, Swing-Away, Power chair covers
Warranty	3-year transferable	3-year transferable