

Options

- **Distance straps** – A variety of Human Care's distance straps in different lengths can be used to achieve a suitable lifting height for both the patient and the caregiver. Must allow minimal distance of 72.8" (1850 mm) between the top position of the bar and the floor
- **Four-point hanger bar** – Allows for versatile and comfortable patient suspension. Can be used with stretcher accessory for horizontal lifting
- **Human Care Data and Service Application** – Collects data and information about the lift's usage to ensure proper use, minimize risk of back injury and increase return on equipment investment. Can access information from PC or tablet device
- **Power traverse** – Allows sideways motorized movement along a rail with a weight capacity of 661 lbs (300 kg), dual speed, high water protection and the possibility of passing curtain gaps in rails
- **Scale** – Medical grade scale can weigh patients up to 705.5 lbs (320 kg)
- **Slings** – A variety of Human Care's slings are compatible with its overhead lift motors to achieve safe, versatile and comfortable patient suspension and transfer. With over 22 different sling models, sized from junior small to bariatric, Human Care slings can be used for:
 - General lifting situations
 - Hygiene lifting in connection to toiling or bathing
 - Repositioning
 - Gait training and rehabilitation
 - Bariatric lifting
 - Specialized lifting situations

Additionally, Human Care slings are available in three different materials, including:

- Polyester
- Polyester net
- Trevira
- **Stretcher** – Light, portable stretcher enables easy horizontal lifting. Must be used with four-point hanger bar.
 - Weight capacity: 441 lbs (200 kg)
 - Unit dimensions: 80.3" x 29.1" x 1.6" (2040 mm x 740 mm x 40 mm)
- **Two-point hanger bar** – Available in four different widths for versatile and comfortable patient suspension:
 - 13.8" (35 cm)
 - 17.7" (45 cm)
 - 21.7" (55 cm)
 - 23.6" (60 cm)

- **Weight capacity** – Three different weight capacities for optimum customization:
 - 331 lbs (150 kg)
 - 485 lbs (220 kg)
 - 661 lbs (300 kg)