

TILT[®]



Toilet Incline Lift

Safety starts at home, and EZ-ACCESS is proud to offer the TILT[®], a patented toilet lift solution that provides dignified bathroom safety and lifts to your natural body movement.



Features & Benefits

- The unit is lightweight, yet strong and durable.
- Installation is easy and typically takes less than 15 minutes.
- Offers great stability, as the TILT attaches directly to the bowl, rather than pushing the unit over the commode.
- Choose between standard or elongated seats.
- Moves the user 7-1/2" forward, which aligns shoulders over feet for optimal balance and positioning.
- Handles remain at the ideal angle, keeping the user's elbow slightly bent for maximum strength positioning.
- The low-profile design allows for convenient lateral and slide-board transfers (in the down position, the arms are lower than the seat).
- The frame is protected by a plastic shield, which can easily be removed for cleaning.
- Unit accommodates users 5' 2" to 6' 4" and bowl heights from 14" to 21".
- Available in corded or battery power.

TILT Toilet Incline Lift



Technical Info

Models:	AC or DC; standard or elongated seat	
Application:	Residential	
Operation:	Unit lifts and lowers users to and from commode	
Mounting:	Bolted directly to bowl	
Cycle time:	Approx. 15 - 20 seconds to lift/lower	
Weight Capacity:	325 lbs.	
Activation:	Wired, push button control with 9-foot coil cord	
Dimensions:	24" W x 18" H x 24" D (approx.)	

	AC (Corded Power)	DC (Battery Power)
Electrical:	10-foot power cord	Rechargeable battery pack
Motor:	120 VAC 60 Hz	24 VDC
Weight:	34 lbs.	40 lbs.

Companion Control

For user and caregiver convenience, the TILT features a hand-held, push button control to operate the lifting process.



Battery Powered Unit

If an outlet is not in close proximity, a battery powered TILT is available.

The rechargeable battery pack mounts to the wall (5" x 16" wall space required) within 3' of the TILT, providing power for approx. 20 cycles with a full charge.

Battery features a 1-year warranty.



Did you know?

According to the CDC, people over the age of 65 suffer the majority of their injuries while getting on and off the toilet.

customerservice@ezaccess.com | 1.800.451.1903 | www.ezaccess.com