

The Jazzy Select® Elite delivers high performance and maneuverability in a compact, stylish design for your active lifestyle. In-line, frontwheel drive technology gives the Jazzy Select Elite excellent stability, provides a tight turning radius around corners for optimal indoor maneuverability, climbs obstacles and grade changes with ease and handles various terrains indoors and out for an overall solid performance. The Jazzy Select Elite comes standard with Pride's exclusive Jazzy® Armor® for today's rental environment,

Jazzy Select® Elite shown in Jazzy Red with high-back seat and Jazzy® Armor® (standard)

depth-adjustable high-back seat, and

flip-up height adjustable foot platform.











SPECIFICATIONS				
<b>Weight Capacity</b>	300 lbs.			
Maximum Speed	Up to 4 mph			
<b>Ground Clearance</b>	1.5"			
Turning Radius	24.75"			
Overall Length	34.75"			
Overall Width	22.75"			
Drive Wheels	9" solid			
Front Wheels	3" solid anti-tips			
Rear Wheels	6" solid casters			
Suspension	Limited			
Drivetrain	Two-motor, in-line, front-wheel drive			
Braking System	Regenerative and electro-mechanical			
Standard Electronics	40 amp, PG GC 3 controller			
Battery Charger	Off-board, 3 amp			
Per-Charge Range	Up to 15 miles			
<b>Battery Requirements</b>	(2) 12 volt, deep cycle, U-1			
Battery Weight	24.5 lbs. each			
Base Weight	68.5 lbs.			
Standard Seat Weight	43 lbs. (comfort, high-back)			
Standard Seat	18" x 18-20"			
Max. Seat Dimensions	20" x 20"			
Warranty	5-year limited on frame; 13-month limited on electronics; 13-month limited on drive motors			

-1	M	u	D	E	L	S
			_	_	_	$\sim$

Jazzy Select Elite-2S-SSCode: K0822Jazzy Select Elite-2S-CCode: K0823

## **FEATURES**

- · Color-through shroud
- Jazzy® Armor® includes shroud and controller quards to protect from daily wear and tear
- · Depth adjustable, black high-back seat
- · Universal, fixed, in-line joystick mount
- Flip-up height adjustable foot platform
- · Black, flat-free, non-marking tires
- 60" lap belt

## **OPTIONS & ACCESSORIES**

- Cup holder
- Weather cover
- Rear basket
- · Oxygen tank holder
- Walker holder
- Cane/crutch holder
- · Saddle bag
- Heel loops

- Lap belt (50", 70", 80", 90")
- Swing-away footrest
- Elevating leg rests
- Military patches
- Swing-away joystick
- · XLR USB charger





**COLOR CHOICES** 

# **OPTIONAL SEATING**

- Synergy<sup>®</sup> Seat
- Solid seat pan

## **RECOMMENDED LIFTS**

- · Backpacker® Plus
- Backpacker® AVP
- · Outlander Full Platform
- Olympian
- Commander



Use the Pride Lift Finder to find a compatible lift for your mobility product:

prideliftsandramps.com



(US) 800-800-8586 • (CAN) 888-570-1113 • pridemobility.com







